Notice d'instruction Instructions for use

Bio
Composants
Médicaux

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Marquage CE : 2013

FIBERFORCE/ FIBER FORCE CST
Instructions for use

Description

FIBERFORCE is a range of glass fibre reinforcement products pre-impregnated with a light-curing methacrylate resin. They may be utilised both in the surgery and in the laboratory, and they provide a solid, visually pleasing and easy-to-use reinforcement. FIBERFORCE is compatible with all types of methacrylate resin (cold curing, hot curing, pressure curing and microwave curing). Preformed items are available in the form of white or pink braid, UD fibre and ribbon in precut sizes. The resin alone is also available in syringes, in white (composite glue) or pink (for repairing acrylic denture plates). FIBERFORCE preformed elements and resins are supplied unpolymerized and are curable in blue light (460-480 nm).

Composition

Preformed glass fibre elements: dimethacrylate urethane (UDMA), triethylene glycol dimethacrylate (TEGDMA), inorganic loading matter, pigments, catalysts/stabilisers.
Resin in a syringe: dimethacrylate urethane (DMAU), triethylene glycol dimethacrylate (TEGDMA), inorganic loading matter, pigments, catalysts/stabilisers.

Indications

- ❖ For use as a structural reinforcement system for manufacture or repair of methacrylate resin dentures.

FIBERFORCE preformed elements should be used in accordance with the application as recommended in the instructions given below. The white or pink resin in a syringe is designed for use with one or more preformed elements.

Contra-indications

Allergy to methacrylates.
Use on a patient in bone growth phase.

Adverse effects

Allergy to methacrylates.

Precautions

- Wear gloves and avoid handling preformed structures with the bare fingers.
- Avoid contact between the unhardened material and the mucosa.
- Unhardened methacrylate resin may cause irritation of the eyes. There may also be a risk of allergy through prolonged contact with this material.
- All work on preformed units should be performed outside the patient's mouth.
- Do not carry out any hardening not described in the instructions (i.e. using the requisite wavelength of blue light) in order to ensure the efficiency of the preformed element.
- When cutting to length or, if necessary, grinding of hardened preformed elements, protective glasses should be worn, together with a breathing mask and gloves, and suitable aspiration should be performed since projected glass fibre particles can irritate the skin. In the event irritation, discomfort may be relieved by use of water and soap (mechanical cleansing).
- Once the packaging has been opened, the pre-impregnated article must be used quickly without exposure to light sources (e.g. sunlight) to avoid all risk of premature hardening.

Use of FIBERFORCE structures

A. Manufacturing methacrylate resin base plates: injection technique

- Anchor the FIBERFORCE mesh firmly in the plaster model to prevent it from moving during injection. This may be achieved as follows:
 - upper jaw: create a small anchor hole in the centre of plaster model (depth: 2 mm) at the level of the papillae and two other holes next to the tuberosities (second molar),
 - lower jaw: create 2 anchor holes in the crests (premolar).
- Apply isolating wax to the plaster model in accordance with standard practice (including in the anchor holes); allow to dry completely.
- Apply a 0.5 mm band of SplintWax to the plaster model and remove excess.
 - Upper jaw: create 6 holes in the wax along the gingival crest (17, 14, 24, 27 and before the Postdam and papillae). These holes will be used to create thickness wedges needed to position the FIBERFORCE mesh properly during injection.
 - Lower jaw: create 5 holes at 37, 34, (31-41), 44 and 47.

- Working away from light, fill the holes created in the wax plaque with FIBERFORCE pink resin: taking care to apply a slight excess of resin.
- Working away from light, take the FIBERFORCE mesh out of its blister strip, cut approximately to size while still in the protective film (using suitable scissors) and then remove the film.
- Isolate the wax with an isolating material for plaster.
- Place the FIBERFORCE mesh in the wax covering the plaster model making sure that it is properly centred in the model.
- Transfer the entire unit to a vacuum forming device and adjust the height of the model in accordance with the instructions for use of the device. Close and create a vacuum to seal the FIBERFORCE mesh on the model. Harden the FIBERFORCE mesh using a photopolymerization chamber, taking care to observe the recommended polymerization times (see end of document).
- Remove the model from the device. Separate the hardened mesh from the model and completely remove all residual wax from the model and the mesh (using boiling water or steam). Allow to dry completely.
- Ensure that the spaces within the mesh are resin free. If necessary, they may be pierced with a handheld tool.
- If necessary, remove any excess by trimming around the mesh: for a more pleasing visual appearance, the mesh should be at least 2 mm from the outer edge. Use a suitable pair of short blunt scissors. Put the mesh to one side.
- Take the plaster model and apply isolating wax in accordance with standard practice; allow to dry completely.
- Prepare the wax using your usual technique. A central injection channel must be prepared in the upper jaw and two lateral channels in the lower jaw. Prepare your muffle in accordance with standard practice.
- Before closing the muffle, fill the anchor holes created in step A1 with FIBERFORCE pink resin, applying an excess of the latter.
- Place the hardened FIBERFORCE mesh on the plaster model in contact with the resin previously included in the anchor holes; the resin must cover the mesh locally. If necessary, place a few drops of resin on top.
- Harden the FIBERFORCE resin in the anchor holes using a photopolymerization chamber: follow the recommended proposition times (see end of document). The mesh is now anchored in position on the plaster model.
- Carry out the injection using the resin you normally use for methacrylate denture plates in accordance with the manufacturer's recommendations.
- Finish the denture in accordance with standard procedure. The FIBERFORCE mesh must be totally encased in the center of the resin.

B. Manufacture of methacrylate resin base plates: pressed technique

- Cover the plaster model with isolating wax in accordance with standard practice and allow to dry completely.
- Apply a 0.5 mm band of SplintWax self-adhesive wax to the plaster and cut around the edges.
 - Upper jaw: create 6 holes in the wax plaque along the gingival crest (17, 14, 24, 27 and before the Postdam and papillae). These holes will be used to create thickness wedges needed to position the FIBERFORCE mesh properly during pressing.
 - Lower jaw: create 5 holes at 37, 34, (31-41), 44 and 47.
- Working away from light, fill the holes created in the wax plaque with FIBERFORCE pink resin taking care to apply a slight excess of resin.
- Working away from light, take the FIBERFORCE mesh out of its blister strip, cut approximately to size while still in the protective film (using suitable scissors) and then remove the film.
- Isolate the wax with an isolating material for plaster.
- Place the FIBERFORCE mesh in the wax covering the plaster model making sure that it is properly centred in the model.
- Transfer the entire unit to a vacuum forming device and adjust the height of the model in accordance with the instructions for use of the device. Close and create a vacuum to seal the FIBERFORCE matrix on the model. Harden the FIBERFORCE mesh in the device using a photopolymerization chamber: taking care to observe the recommended polymerization times (see end of document).
- Remove the model from the device. Separate the hardened mesh from the model and completely remove all residual wax from the model and the mesh (using boiling water or steam). Allow to dry completely.
- Ensure that the spaces within the mesh are resin free. If necessary, they may be pierced with a handheld tool.
- If necessary, remove any excess by trimming around the mesh: for a more pleasing visual appearance, the mesh should be at least 2 mm from the outer edge. Use a suitable pair of short blunt scissors. Put the mesh aside.
- Take the plaster model and apply isolating wax in accordance with standard practice; allow to dry completely.
- Shape your wax and plaster muffle using your normal technique.
- At the time of preparation of the muffle, place the hardened FIBERFORCE mesh on the plaster model. Apply your standard resin for the creation of methacrylate denture plates and press in accordance with the manufacturer's recommendations.
- Finish the denture in accordance with standard procedure. The FIBERFORCE mesh must be totally encased in the resin (horizontally in the centre of the resin).

C. Manufacturing dental splints

- Apply a band of SplintWax Clinique self-adhesive wax to the tracing outlining the position of the future splint.

2. Spread a layer of at least 8 mm of transparent SplintMatrix silicone using the flat mixer tip provided for this purpose. Insert the MatrixTip and cover with more silicone. Allow to harden for the stipulated time (see the instructions for use for the silicone).
3. Remove the silica matrix and discard the excess: it must fit readily in the model. The back of the silicone matrix contains a negative groove, reflecting the SplintWax Clinique wax.
4. Put the SplintWax Clinique wax band to one side.
5. Isolate the model.
6. Working away from light, use suitable scissors to cut out the length of FIBERFORCE braid through the protective film the same length as the wax band.
7. Place a film of composite glue (e.g. white FIBERFORCE resin) in the groove and lay out the braid in the retention groove.
8. Replace on the model without sliding and press.
9. The braid is thus compacted with the dental surfaces to form a ribbon.
10. Extract and discard excess resin removed from the braid using a slowly-turning abrasive paper disc. Sand lightly inside the back.
11. Replace in the silicone matrix.
12. The splint may thus be delivered in its silicone matrix ready for affixation at the dental surgery.

D. Manufacturing dental splints adding one or two teeth (incisor-canine region)

1. Preparation of the artificial tooth

- 1.1 Adjust an acrylic tooth or the composite in the toothless space and lightly sand the palatine or lingual area.
- 1.2 Carry out negative retention in the lingual or palatine aspect where the temporary bridge is to be placed and along 2/3 of its thickness.
- 1.3 Create a counter groove in the direction of the neck.
- 1.4 Applying adhesive (bonding) (e.g. BioBond SF).

2. Preparation of the model

- 2.1 Position the artificial tooth on the model.
- 2.2 Apply a band of SplintWax Clinique self-adhesive wax on the imprint made by the future temporary bridge and push it into the bottom of the horizontal retention groove.
- 2.3 Spread a layer of at least 8 mm of transparent SplintMatrix silicone using the flat mixer tip specially provided. Insert the MatrixTip and cover with more silicone. Allow to harden for the stipulated time (see the instructions for use for the silicone) ensure proper penetration into the retention groove for the artificial tooth.
- 2.4 Remove the silica matrix and discard the excess: it must fit readily in the model. The back of the silicone matrix contains a negative groove, reflecting the SplintWax Clinique wax, with the groove forming a positive excrescence.
- 2.5 Put the SplintWax Clinique wax band to one side.
- 2.6 Isolate the model.
- 2.7 Working away from light, use suitable scissors to cut out the length of FIBERFORCE braid through the protective film the same length as the wax band.
- 2.8 Place a film of composite glue (e.g. white FIBERFORCE resin) in the groove and lay out the braid in the retention groove.
- 2.9 Replace on the model without sliding and press firmly.
- 2.10 Fill the groups in the artificial tooth with FIBERFORCE white resin and apply in the free space.
- 2.11 The braid is thus compacted with the dental surfaces to form a ribbon. It is encased within the horizontal groove.
- 2.12 Fill the palatine depression of the artificial tooth with more composite flow.
- 2.13 Extract and discard excess resin removed from the braid using a slowly-turning abrasive paper disc. Sand lightly inside the back. Replace on the model.
- 2.14 The splint bridge is ready for affixation at the dental surgery.

E. Creating a shaped glass-fibre structure to strengthen an acrylic bridge implant (CST)

Caution: maximum distal cantilever extension on distal implant is 11mm.

1. The plaster laboratory models are fitted with implant analogues and with analogues of conical implant pillars supplied by the implant manufacturer. The latter are used to hold independent posts to be incorporated in the acrylic implant.
2. Isolate the model, install the movable titanium posts: sand, silane and coat with a resin binder (BioBond SF).
3. In the plaster laboratory model create two holes 10 mm from each post at the most distal point: fill the holes with white BioFlow BF light-curing resin and insert a Ø2mm composite fibre traction pin. Light cure.
4. Base braid 1: apply a hybrid 1:6 braid measuring 450 mm to the right tension pin and unwind making one wind anticlockwise around each titanium post, working towards the left tension pin. Maintain under tension and light cure.
5. Base braid 2: wrap around the right tension pin and unwind, making one wind in the same direction around each connection post, and work towards the left tension pin. Maintain under tension and light cure.
6. Transverse braid 3: with the braid still tight, wrap around the right tension pin and unwind, making one anticlockwise wind around each connecting pin, and working towards the left tension pin. Maintain under tension and light cure.

7. Supporting braids 4 and 5: starting from the right tension pin, roll a 1:4 hybrid braid measuring 300 mm around the armature formed by base braids 1 and 2 and transverse braid 3, moving towards the left tension post, then work back in the opposite direction towards the right tension pin. Maintain under tension and light cure.
8. Inject or press the acrylic resin as described in A17 or B13.

F. Repair and strengthening of removable methacrylate resin dentures

1. Inspect the broken denture, and remove any broken pieces if necessary.
2. Glue the denture: clean the fracture line and coat with adhesive (bonding) and FIBERFORCE pink resin. Assemble and photopolymerize.
3. Apply a layer of SplintMatrix transparent silicone at least 8 mm thick inside the back of the broken plate next to the fracture zone. Allow to harden for the stipulated time (see the instructions for use for the silicone). Remove the silicone matrix.
4. Create a transverse groove in the denture at an angle of between 30° and 90° to the fracture line (fissure bur, 1.5 to 2 mm).
5. Replace the transparent silicone matrix.
6. Apply a small quantity of pink FIBERFORCE resin in the groove.
7. Away from light, remove a braid or a UD FIBERFORCE fibre from the blister pack and rapidly trim to the required length using suitable scissors without removing the protective film. When this operation is complete, remove the protective film.
8. Place the braid or UD fibre at the bottom of the groove ensuring that it is well encased in the resin. Several lengths of braid or UD fibre may be used together to ensure optimal reinforcement.
9. Completely fill the groove with FIBERFORCE resin taking care to completely cover the braid or UD fibre. A surplus of resin should be applied.
10. Harden the assembly in a photopolymerization chamber or using a photopolymerization lamp: observe the recommended polymerization times (see below).
11. Polish and finish the assembly in accordance with standard practice. The FIBERFORCE should be totally encased in the resin.

N.B.: You may also use your usual self-curing methacrylate resin: first harden the braid or the UD FIBERFORCE fibre at the bottom of the groove using a polymerization chamber or a polymerization lamp taking care to comply with the recommended polymerization times. Remove it and put aside. Place a little self-curing methacrylate resin at the bottom of the groove. Insert the hardened braid or UD fibre. Then cover the FIBERFORCE reinforcement Splint™ and the self-curing methacrylate resin groove, and harden the resin in accordance with the manufacturer's instructions.

G. Strengthening a denture extension

1. Create a groove in the back of the denture, close to the dentition base, on each side of the extension zone (fissure bur, 1.5 to 2 mm, to a depth of 1 to 2 mm).
2. Place a small quantity of pink FIBERFORCE resin at the bottom of the groove.
3. Away from light, remove a braid or a UD FIBERFORCE fibre from the blister pack and rapidly trim to the required length using suitable scissors without removing the protective film. When this operation is complete, remove the protective film.
4. Place the braid or UD fibre at the bottom of the groove ensuring that it is fully covered by the resin. Prepare a loop in the extension zone: this will be used to do retain the tooth to be added. Several lengths of braid or UD fibre may be used together to ensure optimal reinforcement.
5. Completely fill the groove with FIBERFORCE resin taking care to completely cover the braid or UD fibre. A surplus of resin should be applied.
6. Harden the assembly in a photopolymerization chamber or using a photopolymerization lamp: observe the recommended polymerization times (see below).
7. Cover the loop with FIBERFORCE resin. During the same procedure, affix the extension tooth.
8. Harden the assembly in a photopolymerization chamber or using a photopolymerization lamp: observe the recommended polymerization times (see below).
9. Polish and finish the assembly in accordance with standard practice. The FIBERFORCE reinforcement should be totally encased in the resin.

H. Strengthening of temporary methacrylate resin bridges

1. Remove the temporary methacrylate resin bridge. Create a groove measuring 1 to 2 mm on the occlusal surface along the entire length of the bridge (fissure bur, 1 to 2 mm).
2. Place a little FIBERFORCE white resin in the bottom of the groove.
3. Away from light, remove a braid or a UD FIBERFORCE fibre from the blister pack and rapidly trim to the required length using suitable scissors without removing the protective film. When this operation is complete, remove the protective film.
4. Place the braid or UD fibre at the bottom of the groove ensuring that it is well encased in the resin. Several lengths of braid or UD fibre may be used together to ensure optimal reinforcement.
5. Completely fill the groove with FIBERFORCE resin taking care to completely cover the braid or UD fibre. A surplus of resin should be applied.
6. Harden the assembly in a photopolymerization chamber or using a photopolymerization lamp: observe the recommended polymerization times (see below).
7. Polish and finish the assembly in accordance with standard practice. The FIBERFORCE reinforcement should be totally encased in the resin.

N.B.: You may also use your usual self-curing methacrylate resin: first harden the braid or the UD FIBERFORCE fibre at the bottom of the groove using a polymerization chamber or a polymerization lamp taking care to comply with the recommended polymerization times (see below). Remove it and put aside. Place a little self-curing methacrylate resin at the bottom of the groove. Insert the hardened braid or UD fibre. Then cover the FIBERFORCE reinforcement splint and the self-curing methacrylate resin groove, and harden the resin in accordance with the manufacturer's instructions.

I. Strengthening of wax validation cosmetic dental frame (wax boxing)

1. Apply isolating wax to the plaster model in accordance with standard practice and allow to dry completely.
2. Away from light, remove a white FIBERFORCE braid the blister pack and rapidly trim to the required length using suitable scissors without removing the protective film. When this operation is complete, remove the protective film.
3. On the plaster model, when the braid around the bridge posts, halfway up and crossing over.
4. Harden the assembly in a photopolymerization chamber or using a photopolymerization lamp: observe the recommended polymerization times.
5. Create your wax box in accordance with current recommendations, and include the FIBERFORCE braid.


J. Manufacture of implant-borne removable full dentures: injection technique

1. The mesh must be firmly anchored in order to prevent movement during injection. To do this, create two small anchor holes in the middle of the plaster model (depth: 2 mm).
2. Apply a 0.5 mm band of SplintWax self-adhesive wax to the plaster model. Create 8 holes in the wax plaque along the gingival crest and at the base of the model: these holes will be used to create thickness wedges needed to position the FIBERFORCE mesh properly during injection.
3. Working away from light, fill the holes created in the wax plaque with FIBERFORCE pink resin: taking care to apply a slight excess of resin.
4. Working away from light, take the FIBERFORCE mesh out of its blister strip, cut approximately to size while still in the protective film and then remove the film.
5. Isolate the wax with an isolating material for plaster.
6. Place the FIBERFORCE mesh in the wax covering the plaster model making sure that it is properly centred in the model.
7. Transfer the entire unit to a vacuum forming device and adjust the height of the model in accordance with the instructions for use of the device. Close and create a vacuum to seal the FIBERFORCE mesh on the model. Harden the FIBERFORCE mesh in the device using a photopolymerization chamber, taking care to observe the recommended polymerization times below. It is preferable to use an LED polymerization chamber which does not heat up, and does not melt the wax, thus avoiding the need for laborious cleaning of the mesh after polymerization.
8. Remove the model from the device. Separate the hardened mesh from the model and completely remove all residual wax from the model and the mesh (using boiling water or steam). Allow to dry completely.
9. Check that the holes in the mesh within the mesh are free of resin. If necessary, they may be pierced with a handheld tool.
10. If necessary, remove any excess by trimming poorly adjusted regions around the mesh: for a more pleasing visual appearance, the mesh should be at least 3 mm from the future outer edge. Any pair of short blunt scissors may be used. Follow the edge of the mesh and then place the mesh aside.
11. Take the plaster model and apply isolating wax in accordance with standard practice; allow to dry completely.
12. Prepare the wax using your usual technique. Prepare two injection channels on either side of the anchor holes. Create your muffle in accordance with standard practice.
13. Before closing the muffle, fill the anchor holes created in step J1 with FIBERFORCE pink resin, applying an excess of the latter.
14. Place the hardened FIBERFORCE mesh on the plaster model in contact with the resin previously included in the anchor holes; the resin must cover the mesh locally. If necessary, place a few drops of resin on top.
15. In order to ensure optimal strength of the attachments, apply one or two FIBERFORCE to the mesh, as attachments to the lingual or palatine areas, depending on occlusion conditions.
16. Apply a line of FIBERFORCE pink resin along the braid in order to optimise attachment of the latter to the mesh.
17. Harden the FIBERFORCE braid or braids and resin in the anchor holes using a photopolymerization chamber: observe the recommended polymerization times given below. The mesh is now anchored in position on the plaster model.
18. Carry out the injection using the resin you normally use for methacrylate denture plates in accordance with the manufacturer's recommendations.
19. Polish and finish the assembly in accordance with standard requirements. The FIBERFORCE must be fully immersed in the resin.
20. Affix the attachments either at the chairside or in the laboratory in accordance with the manufacturer's recommendations.

K. Manufacture of implant-borne removable full dentures: pressed technique

1. Using isolating wax, cover the plaster model containing attachment lugs in accordance with standard practice and allow to dry completely.
2. Heat and apply a 0.5 mm wax band to the plaster model. Create 6 holes in the wax plaque along the gingival crest and at the centre of the model: these holes will be used to create thickness wedges needed to position the FIBERFORCE mesh properly during pressing.

3. Working away from light, fill the holes created in the wax plaque with FIBERFORCE pink resin taking care to apply a slight excess of resin.
4. Working away from light, take the FIBERFORCE mesh out of its blister strip, cut approximately to size while still in the protective film and then remove the film.
5. Isolate the wax with an isolating material for plaster.
6. Place the FIBERFORCE mesh in the wax covering the plaster model making sure that it is properly centred in the model.
7. Transfer the entire unit to a vacuum forming device and adjust the height of the model in accordance with the instructions for use of the device. Close and create a vacuum to seal the FIBERFORCE mesh on the model. Harden the FIBERFORCE mesh in the device using a photopolymerization chamber, taking care to observe the recommended polymerization times below.
8. Remove the model from the device. Separate the hardened mesh from the model and completely remove all residual wax from the model and the mesh (using boiling water or steam). Allow to dry completely.
9. Check that the holes in the mesh within the mesh are free of resin. If necessary, they may be pierced with a handheld tool.
10. If necessary, remove any excess by trimming poorly adjusted regions around the mesh: for a more pleasing visual appearance, the mesh should reach the gingival border but without covering the vestibular aspect. Any pair of short blunt scissors may be used. Return the polymerized mesh to the plaster model
11. In order to ensure optimal strength of the attachments, apply one or two FIBERFORCE to the mesh, as attachments to the lingual or palatine areas, depending on occlusion conditions.
12. Apply a line of FIBERFORCE pink resin along the braid in order to optimise attachment of the latter to the mesh.
13. Harden the FIBERFORCE braid or braids and resin in the anchor holes using a photopolymerization chamber: observe the recommended polymerization times given below. Put the mesh aside.
14. Take the plaster model and apply isolating wax in accordance with standard practice; allow to dry completely.
15. Create your wax and muffle using your usual techniques.
16. When closing the muffle, place the hardened FIBERFORCE mesh on the plaster model. Cover with the resin you normally use for methacrylate resin dentures and press in accordance with the manufacturer's recommendations.
17. Polish and finish the assembly in accordance with standard requirements. The FIBERFORCE mesh should be totally encased in the resin.
18. Affix the attachments either at the chairside or in the laboratory in accordance with the manufacturer's recommendations.

Store at between 12 and 30°C in closed packaging.

Can cause skin allergies.

Date of last leaflet update: 03.2020

CE marking: 2013

Polymerization times for FIBERFORCE products

In the surgery, using a hand-held photopolymerization lamp:

Lamp type	Blue LED, 5W	Halogen, 1100 mW/cm ²	Halogen, 550 mW/cm ²
Minimum time	20 s	30 s	80 s

In the laboratory, in a photopolymerization chamber:

Lamp type	Stroboscopic xenon, 250 mW/cm ²	Neon, 6800 mW/cm ²	Mercury vapour (arc)
Minimum time	4 min	10 min	20 min

Lamp type	Halogen, 550 mW/cm ²	Halogen, 1100 mW/cm ²	Blue LED
Minimum time	2 min	40 s	5 to 10 min (depending on power)

N.B.:

We strongly advise against use of PEI chambers (fitted for example with 4 x 9W Radium Ralutec or Osram Dulux S blue tubes) due to incompatibility with the camphorquinone absorption spectrum.

Note

This product is designed solely for use in dental applications.
Keep out of the reach of children.

Once polymerized, the device cannot be reused.

Do not reuse the tips for syringes as there is a risk of cross contamination.

Single-use product.

Store protected from light.

NB: see instructions for use.