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Réf. : 1 04 001 034 B

FIBERFORCE D-Lab
Instructions for use

Description

Fiber Force D-Lab is an accurate system to relate implant position using a high stability light-curing preimpregnated fiber. For chairside procedures, preformed glass fiber elements of precut lengths that will be incorporated into a translucent dental silicone impression. Their blue color makes them easy to locate through the translucent silicone. The unfilled resin is available in refill syringes. FIBER FORCE preformed elements and resins are supplied unpolymerized and are curable using blue light (460-480 nm).

Composition

Preformed glass fibre elements: glass fiber, urethane dimethacrylate (UDMA), triethylene glycol dimethacrylate (TEGDMA), inorganic fillers, pigments, catalysts/stabilizers.

Resin syringe: urethane dimethacrylate (UDMA), triethylene glycol dimethacrylate (TEGDMA), inorganic fillers, pigments, catalysts/ stabilizers.

Indications

Positioning system for transparent silicone in implant transfer impressions.

Contra-indications

Allergy to methacrylates.

Adverse effects

To date, there are no known adverse effects.

Precautions

- Avoid any contact of the oral mucosa with uncured material.
- Wear gloves and avoid handling pre-cut fibres with bare hands.
- The protocol must be respected to guarantee complete polymerization of the precut fibres and to ensure dimensional stability of the system.
- Uncured methacrylate resin may cause irritation of the eyes. There may also be a risk of allergy through prolonged contact with this material.
- For length adjustment before polymerization, please use suitable scissors.
- In the case of cutting polymerized fibres in the mouth, work under irrigation and with surgical aspiration. In the event of irritation, discomfort may be relieved by rinsing with water. Once the packaging has been opened, the pre-impregnated fibres must be used quickly without exposure to light sources to avoid risk of premature curing.
- Cut the first centimeters of the braid before use.

Use of FIBERFORCE elements

Prior to use:

1. The implant transfer impression copings are placed by the practitioner.
2. The transfer copings are isolated and the fixation screws are protected by plastic covers or other isolators.
3. Cut the first centimeters of the braid.

D-Lab system

4. A 750 mm 1:4 hybrid braid is clamped distally to the most distal right-hand transfer coping and a wrapped twice around the coping.

5. Base braid 1: unwind the hybrid braid counter clockwise wrapping around each transfer coping, working towards the more distal left coping. Maintain tension and light cure. Maintain tension by pressing the button and light cure.

6. Base braid 2: Wrap around the more distal left coping and unwind, wrapping around each coping towards the right transfer coping. Maintain tension by pressing the button and light cure.

7. Upper braid 3: With the braid still under tension, unwind the braid counter clockwise, wrapping once around each transfer coping, working towards the more distal left coping. Make a final turn around this last coping. Maintain tension by pressing the button and light cure. Then cut the fiber with scissors.

8. Apply light-cure resin to the areas where fibres meet the copings and light cure. Carefully light cure the whole assembly for 40 seconds, including the fibre links between the implants.

9. Fill a suitable mesh impression tray with a transparent dental silicone impression material.

10. Insert the impression tray to the prepared zone and maintain pressure while the material is setting.

11. continue to light cure the Fiber Force braid through the translucent silicone.

12. Remove the plastic protections.

13. Unscrew the implant transfer screws and remove the impression tray from the mouth.

14. An immediate load implant prosthesis can now be completed with the CST procedure in a prosthetic laboratory.

Note

For dental use only

Keep out of reach of children.


Single use product

Store protected from light

Caution: See instruction for use

Store between 12 and 30°C, in closed packaging.

Can cause skin allergies.

Eliminate waste according to regulations.

Date of last leaflet update: 09.2020

CE marking : 2019